

| ID | Requirement | Related Use Case | Fulfilled By | Test | Description |
|----|--|------------------|-----------------------------------|--|--|
| 1 | A light on the user interface that indicates an active pulse reading | N/A | mainwindow | Run project and observe the user interface. | A QLabel is used to display different states. |
| 2 | The user interface consists of a screen and buttons. the screen contains the menu and display graph. | N/A | mainwindow | Run project and observe the user interface. | The screen in the mainwindow is a custom widget, which can display menus and graph, and the mainwindow also contains 8 buttons.a menu button, a standard back button which will return the user to the menu, four arrow buttons and a selector in the center of the arrow buttons which also functions as a start/stop button. |
| 3 | A light on the user interface that indicates coherence level | N/A | mainwindow | Start new session and observe the coherence light. | Qpushbutton is used to indicates coherence level, when indicates coherence level changed, the background color of the button will change. |
| 4 | The menu consists of the following options: start new session, settings, log/history | N/A | menu | Run project, press power button, and observe the menu | The menu class contains a QStringList which contains the current menu items, a parent pointer holds the parent Menu and a Qvector holds the sub-Menus. |
| 5 | session screen displays the main HRV graph. | N/A | mainwindow, qcustomplot, hrvgraph | Start new session and observe the screen. | The qcustomplot is used to draw graph. The hrvgraph contains a Qvector contains the graph data. |
| 6 | The settings menu includes challenge level and breath pacer settings. | N/A | menu | press power button and select the settings menu. | The menu class contains a QStringList which contains the current menu items, |
| 7 | There are 4 challenge levels for coherence | N/A | menu | press power button and select the settings menu, then select challenge level menu. | The menu class contains a QStringList which contains the current menu items, |

| | | | | | |
|----|--|-----|--------------------------------|--|---|
| 8 | There are 4 breath pacer frequencies in menu | N/A | menu | press power button and select the settings menu, then select breath pacer settings menu. | The menu class contains a QStringList which contains the current menu items, |
| 9 | The menu contains a log tab of all sessions, when selected show the summary view. | N/A | mainwindow,menu,record,session | press power button and select the log menu. | The mainwindow class contains a QVector which contains all sessions, each session is stored in the form of record. |
| 10 | The user can clear log of all sessions. | N/A | mainwindow,menu,record | select the log menu, then select the clear menu. | when clear menu is selected, the session container in the mainwindow will be clean. |
| 11 | The user can not start new session when HR sensor disconnected. | N/A | mainwindow | disconnect the HR sensor, then select start new session menu. | when start new session menu is selected, the mainwindow will check the value of HrContactComboBox, if value is 0, then display a QMessageBox to remind the user that the HR sensor is disconnected. |
| 12 | A strip of lights shows the breath pacer status. | N/A | mainwindow | Run project and observe the user interface. | A Qprogressbar is used to show the breath pacer status. |
| 13 | A strip of lights shows the power level. | N/A | mainwindow | Run project and observe the user interface. | A Qprogressbar is used to show the battery level. |
| 14 | When the battery is low, a reminder box will pop up. | N/A | mainwindow | start new session, then wait until the battery level is below 10 | when battery level changed, mainwindow will check the value of battery level, if value is below 10, a QMessageBox is pop up to remind user that the battery level is low. |
| 15 | When the HR sensor disconnected in a session, the session will close and a reminder box will pop up. | N/A | mainwindow | start new session, then disconnect the HR sensor. | when sensor connection state changed, the mainwindow will check the value of HrContactComboBox, if value is 0, then display a QMessageBox to remind the user that the HR sensor is disconnected. |

| | | | | | |
|----|---|-----|------------|---|--|
| 16 | When the power is 0, the device cannot be turned on | N/A | mainwindow | Adjust the battery level to 0, then press the power button. | When power button is pressed, mainwindow will check the battery level, if battery level is 0, the device powerStatus will be set to false. |
| 17 | The application does not contain any memory leaks. | N/A | N/A | use mtrace or valgrind to trace the memory allocation. | All dynamically allocated memory that the program was designed to allocate is deleted in the appropriate class destructor. |